



Dr Trust®

USFDA & CE APPROVED



## Dr Trust Advanced Blood Pressure Monitor

Manage The Performance Of Your Heart



*Rohit Sharma*  
Rohit Sharma

Compatible With:  
**Dr Trust 36**



## Dr Trust 36 One Application For All Your Health Needs

- Smart Sensei Technology
- Share Reports With Doctors
- Get Health Recommendations
- Connect All Dr Trust Devices

## Unique Features



Fuzzy Algorithm Makes Monitoring Comfortable



Works With Mobile Charger & Power Bank

## Why Choose Dr Trust?



High Definition Screen



WHO Indicator



Automatic Power Off Feature



Dual Users Mode



Adjustable Cuff



## Contents In The Box

Dr Trust Smart BP Monitor



Four AAA Batteries



Comfortable Cuff (22 x 42cm)



Instructions Manual Card



## Specifications

Measurement Procedure:	Oscillometric, corresponding to Korotkoff method: Phase I: systolic, Phase V: diastolic
Display:	Digital display
Measuring range:	Pressure: 30 to 280 mmHg (in 1 mmHg increment) Pulse: 40 to 199 beat/minute
Static accuracy:	Pressure: ±3mmHg / Pulse: ±5% of reading
Measuring Resolution:	1mmHg
Inflation:	Automatic inflation by internal pump
Memory Function:	2 x 120 memories for 2 users (SYS, DIA, Pulse)
Decompression:	Constant exhaust valve system
Power Source:	4- size "AA" alkaline Batteries (micro USB DC5.0V / 1.0A) optional
Operation Temperature:	5~40°C / 41~104°F
Operation Humidity:	15%~85%RH maximum
Storage Temperature:	-20~60°C / 14~131°F
Storage Humidity:	10%~95%RH maximum
Dimensions:	142 x 98 x 53 ±1.0 mm
Weight:	600 g±5g (including batteries and cuff)
Cuff Pressure Display Range:	0~290mmHg/0~38.7kPa
Electrical Shock Protection:	Internal power unit
Safety Classifications:	Type B equipment
Mode of operation:	Continuous operation
Protection against ingress of water:	IP22
Accessories:	M/L-size Cuff, 4 "AA" batteries, instruction manual